**Warning** :The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Implant Systems. Guidance in the handling of the Dyna Helix® ART Octalock® Implant System and Dyna (Octalock®) Implant System by an experienced operator is strongly recommended. Dyna Helix® ART Octalock® and Dyna (Octalock®) Implant Systems must only be used by properly trained dentists/doctors and in combination with original components. In case of multiple use the following dangers could occur: cross infection, damaging products and as a result their function, wrong identification of products. For more detailed information please refer to the Dyna Implant Manuals as well as Dyna Terms of Guarantee – available on request. With the publication of this instructions for use all previous are no longer valid.

#### Product description

The Dyna Helix® ART Octalock® Implant System consist of titanium cylinder-like, screw type implants with ART (Acid etched Rough Titanium) coating and a range of prosthetic abutments. The implants are available in various lengths and diameters. The Dyna Helix® ART Octalock® Implants can be used in a one or two stage surgical procedure taken the bone quality and the operation technique in consideration. The Dyna (Octalock®) Implant System includes titanium cylinder-like, push-in type implants with HA coating and a range of prosthetic abutments. The implants are available in various lengths and diameters. The systems besides implants contains also surgical, prosthetic and laboratory components.

#### Materials

The following materials Dyna uses for its implant products: Titanium and Hydroxyapatite.

#### Content package

See label on packaging.

#### Indications

The Dyna (Helix® ART) (Octalock®) Implant System is intended for restoring edentulous or partially edentulous patients with fixed, removable or fixed detachable protheses. Dyna Helix® TM implants are specially designed for single-phase technology. The Dyna (Helix ART) Octalock® Implant System is intended for all types of oral implant supported restorations including single missing teeth. The Dyna Helix® ART Octalock® Implant System has been provided with the Dyna Octalock® connection, therefore, all prosthetic components are interchangeable with all Dyna Octalock® Implant Systems.

#### Contraindications

All contraindications associated with elective oral surgery should be noticed. These include: local (e.g. local infection, inadequate hygiene, unfavourable bite relation, mucosa infections, jaw defects, bruxism, periodontal disease etc), systemic (e.g. bone vascular disorders, cardiac disease, nephritis or nephroses, rheumatic disease, defective immune response, allergic diseases etc), age related limitations (e.g. children) and other (e.g. pregnancy, chemotherapy, lack of motivation, diseases treated with large doses of steroids etc.). The use of the Dyna Push-in Type Implant Systems (Dyna Implant System and Dyna Octalock® Implant System) is contraindicated for fixed constructions in combinations with Ø3,0mm implants; immediate implant placement; direct/immediate loading and restoring single missing teeth. The use of the Dyna (Helix® ART) (Octalock®) Implant System is contraindicated, in all situations where it is not indicated.

#### Surgical method (shortened version)

All Dyna implants

► Make sure that the indications for implant placement are met.

► According to clinical situation anaesthetize, make incision and raise the mucosal flap.

► Under excessive cooling (sterile saline solution) prepare the implant site; depth markings (8, 10, 11,5, 13, 15, 17 mm) are clearly indicated on each drill.

► Fixed sequence of drilling is obligatory (Helix®: pilot, ø3,2; ø3,6, ø4,2, ø5,0 mm) (Dyna (Octalock®): pilot, ø3,0; ø3,6, ø4,0, ø5,0 mm) the recommended speed for the implant site preparation with Dyna twist drills is max. 800 rpm (decrease the drilling speed when using larger diameter drills: ø3.6-600rpm, ø4.2-500rpm, ø5.0-400rpm). Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the minimum bone width for ø3.6 implants is 5.6mm).

Dyna Helix®

► Finalize the shape of the site with the Dyna Helix® Cortical Reamer (This does not apply to the Dyna Helix® TM and Dyna Helix® ST implants). Depth according to the bone quality. (Type 1= 6mm, Type 2=4-6mm, Type 3=2-4mm, Type 4=0-2mm) For details see Dyna Helix® ART Octalock® Implant manual.

► According to the bone quality use the Dyna Helix® Taper. In type 1 bone the Taper angles should be used. For details see Dyna Helix® ART Octalock® Implant manual.

All Dyna implants

► Check the site with the Dyna Parallel/depth Instruments/Gauges and perform radiographic examination.

Dyna (Octalock®).

► Open the double peel pouch packaging and remove the implant with the closing screw by gently pulling the top of the tube, insert the implant gently into the preparation and tighten the closing screw.

Dyna Helix®

► Open the double peel pouch packaging and remove the top of the tube (closing screw included). Push the Dyna Octa driver for the Torque Wrench or Latch head into the internal Octagon of the implant and remove the implant gently. Although the implant will touch the titanium spring it will not damage the implant. Insert the implant gently into the preparation; thread it into its final position with maximum 30 rpm and an optimum torque range between 25 and 55Ncm, but never exceeding 60Ncm. Tighten the closing screw with 10Ncm or when indicated the chosen abutment with 30Ncm (note: extension abutments with 35Ncm).

► Prevent (over)loading Dyna Helix® (TM) implants for minimally 8 weeks. When biomechanical criteria for immediate loading are not met this period of time must be extended.

All Dyna implants

► Suture the site. Give the patient indications. If required and necessary make a temporary. In order to check the implants' position and control its healing process it is recommended to perform radiographic examination during and/or after the surgery, and later during the healing period. Implants must achieve sufficient primary stability.

#### Precautions

Proper planning is essential for successful implant treatment.

Beware of undercut planning that the marks on the drills for the Helix® implants are 0,25 mm higher and for the Dyna (Octalock®) push-in implants 0,75mm higher as the indicated implant length.

Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.

Avoid thermal trauma to the bone; use the largest possible diameter of implants.

Use only in combination with original instruments.

Implant mobility, bone loss, infection may be the symptoms of implant failure.

Helix® Implant should be inserted with the optimum torque ranging between 25-55Ncm.

Mechanical damage of the Helix® Implant can occur above a torque of 70Ncm and of latch-head instruments above a torque of 50Ncm

Avoid overheating during insertion of the implant.

Dyna implants must not be altered in any way.

Rotary instruments must not be used more then 20 times, or when they become damaged or dull.

In spite of special treatments for wear resistance instruments can wear after time. Change the instruments for new when there any doubt regarding the fit to prevent damaging implants, abutments, screws, etc.

In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

#### Sterilization

All Dyna implants are sterile, supplied in a double peel-pouch packaging, and intended for single use only. Other Dyna products are supplied non-sterile and must be sterilized or disinfected before clinical use, in an appropriate manner. Under no circumstances sterilize or re-sterilize implants. Do not use implants when the packaging is damaged.

#### Complications

Apart from the general complications associated with any oral surgery the following, related to dental implants, are the most frequent: inflammation, infection, bone loss, swelling, (chronic) pain, paresthesia, patient discomfort, tissue degeneration, bone/implant/restoration fracture, implant mobility, exfoliation, injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Implant Manuals.

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#### Copyright y marcas comerciales