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<p>Product description</p> <p>The Dyna EFM Alloy and Dyna EFM keeper are parts of the Dyna (Direct) System, a magnetic retention system. Dyna EFM alloy is a ferromagnetic alloy to be casted as root keeper for magnets whereas the Dyna EFM keeper is a pre-fabricated, ferromagnetic alloy keeper to be used directly on endodontically treated roots.</p>																																			
<p>Materials</p> <p>Dyna EFM Alloy, also used for fabrication of the pre-fabricated Dyna EFM Keeper, is a white precious alloy consisting of Pd 60,0%, Co 38,0%, Ga 1,2%, Pt 1,0%</p> <p>Further specifications:</p> <table> <tbody><tr> <td>Specific gravity</td> <td>: 10,0 g/cm3</td> <td>0,2% yield strength</td> <td>: 500 N/mm2</td></tr> <tr> <td>Elongation</td> <td>: 15%</td> <td>Total investment expansion</td> <td>: 1,9%</td></tr> <tr> <td>Preheating temperature</td> <td>: 850°C (1562°F)</td> <td>Melting range</td> <td>: 1195-1210°C</td></tr> <tr> <td>Casting temperature</td> <td>: 1350°C (2462°F)</td> <td>Hardness</td> <td>: 295 kg/mm2</td></tr> <tr> <td>Tensile strength</td> <td>: 80 daN/mm2</td> <td>Specific gravity</td> <td>: 10,00 9/cm3</td></tr> <tr> <td>Solder</td> <td>: 1060° - 1070°C</td> <td></td> <td></td></tr> </tbody></table> <p>Magnetic data</p> <table> <tbody><tr> <td>Magnetic permeability</td> <td>: UM 240</td> <td>Flux density</td> <td>: B 10,0 kg</td></tr> <tr> <td>Coercitive force</td> <td>: HC 5,7 Oersted</td> <td></td> <td></td></tr> </tbody></table>	Specific gravity	: 10,0 g/cm3	0,2% yield strength	: 500 N/mm2	Elongation	: 15%	Total investment expansion	: 1,9%	Preheating temperature	: 850°C (1562°F)	Melting range	: 1195-1210°C	Casting temperature	: 1350°C (2462°F)	Hardness	: 295 kg/mm2	Tensile strength	: 80 daN/mm2	Specific gravity	: 10,00 9/cm3	Solder	: 1060° - 1070°C			Magnetic permeability	: UM 240	Flux density	: B 10,0 kg	Coercitive force	: HC 5,7 Oersted					
Specific gravity	: 10,0 g/cm3	0,2% yield strength	: 500 N/mm2																																
Elongation	: 15%	Total investment expansion	: 1,9%																																
Preheating temperature	: 850°C (1562°F)	Melting range	: 1195-1210°C																																
Casting temperature	: 1350°C (2462°F)	Hardness	: 295 kg/mm2																																
Tensile strength	: 80 daN/mm2	Specific gravity	: 10,00 9/cm3																																
Solder	: 1060° - 1070°C																																		
Magnetic permeability	: UM 240	Flux density	: B 10,0 kg																																
Coercitive force	: HC 5,7 Oersted																																		
<p>Content package</p> <p>See label on packaging.</p>																																			
<p>Indications</p> <p>Dyna magnetic systems are intended for increasing retention of (partial) overdentures, supported on implants and natural roots. Magnets are specially indicated for overdentures in patients with manual dexterity deficiency and/or severe reduced alveolar bone.</p>																																			
<p>Contraindications</p> <p>All contraindications associated with elective prosthetic rehabilitation on endodontically compromised teeth should be noticed. The use of Dyna Magnetic systems is additionally contraindicated in situations where:</p> <ul style="list-style-type: none">▶ proper articulation can not be restored ▶ proper alignment of the magnet and keeper can not be guaranteed ▶ patients are known to be allergic to one or several of the elements contained in the attachment materials. In case of doubt perform preliminary allergological testing <p>The use of Dyna magnetic systems is contraindicated in situations where it is not indicated.</p>																																			
<p>Laboratory procedure</p> <p>The Dyna EFM alloy is castable with every casting machine or gas/oxygen torch. Use a clay crucible (non carbon) and flux. Oxidation layer remains on molten alloy. Alloy does not shine at casting temperature. Remove around 0,1mm of upper layer of the keeper while finishing.</p> <p>Recasting : use at least 30% of new material</p> <p>Soldering : as all palladium alloys, see manual for details</p> <p>Each product has its special clinical and laboratory procedure. See Dyna Magnet Manual for detailed information.</p>																																			
<p>Clinical application</p> <p>Whenever working with overdenture attachments take the following into consideration:</p> <ul style="list-style-type: none">▶ allowing the patient 2-4 weeks of wear-in time ▶ type of implant/root-overdenture connection (resilient vs. non-resilient) ▶ alignment of the male and female parts ▶ use Dyna Magnet Bond to realize chemical bond between attachment and acrylic. <p>For details see the Dyna Magnet Manual.</p>																																			
<p>Disinfection</p> <p>Dyna EFM alloy is a dental alloy. Follow the usual pattern of disinfection/sterilization for dental alloys and castings. Dyna EFM keepers are supplied not sterile in a blister packaging, and intended for single use only. They must be disinfected before clinical use in an appropriate manner. Do not use when the packaging is damaged.</p>																																			
<p>Precautions</p> <p>When damaged exchange products immediately.</p> <p>Improper technique can contribute to treatment failure and/or bone loss. Check the overdenture/attachments regularly (every 6 months).</p> <p>Use only the original instruments.</p> <p>Denture cleaning agents should be used with care. Overnight cleaning solutions should be avoided.</p> <p>Magnets must not be heated above 120°C / 248°F and may only be used in combination with Dyna EFM alloy, Dyna Keepers or Dyna Medical (magnet) abutments.</p> <p>In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.</p>																																			
<p>Handling and Storage</p> <p>Store in clean, dry, dust-free, dark room at room temperature.</p>																																			
<p>Delivery</p> <p>Federal law restricts these devices to sale by or on the order of a dentist or a physician.</p>																																			
<p>Traceability of serial/lot numbers</p> <p>It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label does not contain information for traceability purposes, but only the La-beler Identification Code (LIC) "EDYN" followed by the Dyna reference number. Read the human-readable interpretation of the bar code symbol printed below for verification.</p>																																			
<p>Training</p> <p>Dyna Dental Engineering BV arranges regular training courses for the beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of Dyna Implant System.</p>																																			
<p>Copyright and trademarks</p> <p>All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® and Dyna Helix® are registe-red trademarks of Dyna Dental Engineering BV.</p>																																			

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