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#### Schulung

**Warning** :The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Implant Systems. Guidance in the handling of the Dyna Helix® Implant System by an experienced operator is strongly recommended. Dyna Helix® Implant System must only be used by properly trained dentists/doctors and in combination with original components. In case of multiple use the following dangers could occur: cross infection, damaging products and as a result their function, wrong identification of products. For more detailed information please refer to the Dyna Implant Manuals as well as Dyna Terms of Guarantee – available on request.

With the publication of this instructions for use all previous are no longer valid.

**Product description**

The Dyna Helix® Implant System consist of titanium cylinder-like, screw type implants with ART (Acid etched Rough Titanium) coating and a range of prosthetic abutments. The implants are available in various lengths and diameters. The Dyna Helix® Implants can be used in a one or two stage surgical procedure taken the bone quality and the operation technique in consideration. The systems besides implants contains also surgical, prosthetic and laboratory components.

**Materials**

The following materials Dyna uses for its implant products: Titanium.

**Content package**

See label on packaging.

**Indications**

The Dyna Helix® Implant System is intended for restoring edentulous or partially edentulous patients with fixed, removable or fixed detachable prostheses. Dyna Helix® TM implants are specially designed for single-stage technology. The Dyna Helix® Implant System is intended for all types of oral implant supported restorations including single missing teeth. The Dyna Helix® Implant System has been provided with the Dyna Octalock® connection, therefore, all prosthetic components are interchangeable with all Dyna Octalock® Implant Systems.

**Contraindications**

All contraindications associated with elective oral surgery should be noticed. These include: local (e.g. local infection, inadequate hygiene, unfavourable bite relation, mucosa infections, jaw defects, bruxism, periodontal disease etc), systemic (e.g. bone vascular disorders, cardiac disease, nephritises or nephroses, rheumatic disease, defective immune response, allergic diseases etc), age related limitations (e.g. children) and other (e.g. pregnancy, chemotherapy, lack of motivation, diseases treated with large doses of steroids etc.). The use of the Dyna Helix® Implant System is contraindicated, in all situations where it is not indicated.

**Surgical method (shortened version)**

► Make sure that the indications for implant placement are met.

► According to clinical situation anaesthetize, make incision and raise the mucosal flap.

► Under excessive cooling (sterile saline solution) prepare the implant site; depth markings (8, 10, 11.5, 13, 15, 17 mm) are clearly indicated on each drill.

► Fixed sequence of drilling is obligatory (Helix®: pilot, ø3.2; ø3.6, ø4.2, ø5.0 mm). The recommended speed for the implant site preparation with Dyna twist drills is max. 800 rpm (decrease the drilling speed when using larger diameter drills: ø3.6-600rpm, ø4.2-500rpm, ø5.0-400rpm). Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the mini-mum bone width for ø3.6 implants is 5.6mm).

► Finalize the shape of the site with the Dyna Helix® Cortical Reamer (This does not apply to the Dyna Helix® TM and Dyna Helix® ST implants). Depth according to the bone quality. (Type 1= 6mm, Type 2=4-6mm, Type 3=2-4mm, Type 4=0-2mm) For details see Dyna Helix® Implant manual.

► According to the bone quality use the Dyna Helix® Taper. In type I bone the Taper always should be used. For details see Dyna Helix® Implant manual.

► Check the site with the Dyna Parallel/depth Instruments/Gauges and perform radiographic examination.

► Open the implant packaging. The implant is freely seated in the titanium tube (closing screw included). Push the Dyna Octa driver for the Torque Wrench or Latch head into the internal Octagon of the implant and remove the implant gently. Insert the implant gently into the preparation; thread it into its final position with maximum 30 rpm and an optimum torque range between 25 and 55Ncm, but never exceeding 60Ncm. Tighten the closing screw with 10Ncm or when indicated the chosen abutment with 30Ncm (note: extension abutments with 35Ncm).

► Prevent (over)loading Dyna Helix® (TM) implants for minimally 8 weeks. When biomechanical criteria for immediate loading are not met this period of time must be extended.

► Suture the site. Give the patient indications. If required and necessary make a temporary. In order to check the implants' position and control its healing process it is recommended to perform radiographic examination during and/or after the surgery, and later during the healing period. Implants must achieve sufficient primary stability.

**Precautions**

Proper planning is essential for successful implant treatment.

Beware during implant planning that the marks on the drills for the Helix® implants are 0.25 mm higher as the indicated implant length.

Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care. Avoid thermal trauma to the bone; use the largest possible diameter of implants.

Use only in combination with original instruments.

Implant mobility, bone loss, infection may be the symptoms of implant failure.

Helix® Implant should be inserted with the optimum torque ranging between 25-55Ncm.

Mechanical damage of the Helix® Implant can occur above a torque of 70Ncm and of latch-head instruments above a torque of 50Ncm

Avoid overheating during insertion of the implant.

Dyna implants must not be altered in any way.

Rotary instruments must not be used more than 20 times, or when they become damaged or dull.

In spite of special treatments for wear resistance instruments can wear after time. Change the instruments for new when there any doubt regarding the fit to prevent damaging implants, abutments, screws, etc.

In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

**Sterilization**

All Dyna implants are sterile, supplied in a double peel-pouch packaging, and intended for single use only. Other Dyna products are supplied non-sterile and must be sterilized or disinfected before clinical use, in an appropriate manner. Under no circumstances sterilize or re-sterilize implants. Do not use implants when the packaging is damaged.

**Complications**

Apart from the general complications associated with any oral surgery the following, related to dental implants, are the most frequent: inflammation, infection, bone loss, swelling, (chronic) pain, paresthesia, patient discomfort, tissue degeneration, bone/implant/restoration fracture, implant mobility, exfoliation, injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Implant Manuals.

**Please note:**

It is the user of Dyna products who is obliged to determine whether or not any products are suitable for a particular clinical situation. It is the user of Dyna products who is obliged to document in appropriate manner the products used for each patient. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications. When using our product intra-orally take proper care to prevent them from being inhaled or ingested.

**Handling and Storage**

Store in clean, dry, dust-free, dark room at room temperature. Do not use after expiry date indicated on the packaging.

**Delivery**

Federal law restricts these devices to sale by or on the order of a dentist or a physician.

**Traceability of serial/lot numbers**

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label does not contain information for traceability purposes, but only the Labeler Identification Code (LIC) "EDYN" followed by the Dyna reference number. Read the human-readable interpretation of the bar code symbol printed below for verification.

**Training**

Dyna Dental Engineering BV arranges regular training courses for the beginning and advanced implanto-logists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of Dyna Implant System.

**Copyright and trademarks**

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authori-zation of Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV..

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