ENGLISH - Instructions for use DYNA ATTACHMENTS

Warning

The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Systems. Guidance in the handling of the Dyna (Direct) System, Dyna Helix® Implant System by an experienced operator is strongly recommended. Dyna Helix® Implant Systems must only be used by properly trained dentists/doctors and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the Dyna Manuals as well as Dyna Terms of Garantie – available on request. With the publication of this instructions for use all previous versions are no longer valid.

Product description

Dyna attachments include the Magnetic attachment system – magnetically retained prostheses

Materials

Dyna uses for its attachments: case of Surgical Steel (316L) and magnet 26%Nd-70%Fe-1%B

Content package

See label on packaging.

Indications

Dyna attachments are intended for increasing retention of (partial) overdentures, supported on implants and natural roots in combination with magnet abutments. Magnets are specially indicated for overdentures in patients with manual dexterity deficiency and/or severely reduced alveolar bone.

Contraindications

All contraindications associated with elective oral surgery and prosthetic rehabilitation on implants should be noticed. The use of Dyna attachments (in combination with magnet abutments) is additionally contraindicated in situations where:

- ▶ proper articulation can not be restored
- ▶ proper alignment of the male and female parts can not be guaranteed
- ▶ patients known to be allergic to one or several of the elements contained in the attachment materials.

In case of doubt perform preliminary allergological testing

The use of Dyna attachments is contraindicated in situations where it is not indicated.

Laboratory procedure

Each attachment has its own special fixation procedure. See the Dyna manuals for detailed information.

Clinical application

Whenever working with overdenture attachments take the following into consideration:

- ▶ allowing the patient 2-4 weeks of wear-in time
- ▶ type of implant/root-overdenture connection (resilient vs. non-resilient)

For details see the Dyna Manuals.

Disinfection

Dyna attachments are supplied not sterile in a blister packaging, and intended for single use only. They must be disinfected in an appropriate manner before clinical use. Do not use when the packaging is damaged.

Precautions

When damaged exchange attachments immediately. Do not grind or sandblast. Improper technique can contribute to implant/prosthesis failure and/or bone loss. Check the overdenture/ attachments regularly (every 6 months). Use only the original instruments.

Denture cleaning agents should be used with care. Overnight cleaning solutions should be avoided. Implant mobility, bone loss, infection may be the symptoms of implant failure. Magnets may not be heated above 120°C / 248°F and only be used in combination with Dyna EFM alloy, Dyna Keepers and Dyna Medical(magnet) abutments. Denture without balanced occlusion and articulation will contribute to failure.

When treating an implant patient, facilitate a disinfected or sterile environment.

Facilitate safe disposal of devices according to safety precautions in the general dental practice or hospitals.

In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

Magnetic fields

The research on influence of the magnetic fields on mammalian tissues shows that using permanent mini-magnets for dental purposes poses no threat to the patients. Magnets may influence the image of MRI especially when it concerns the head and neck region. Magnets and magnetic abutment keepers should be removed before the scan. The influence of the magnets on pacemakers devices is limited and when magnets are used according to their indications negligible.

Please note:

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority. When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested.

Handling and Storage

Store in clean, dry, dust-free, dark room at room temperature.

Delivery

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.**Traceability of serial/lot numbers**

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) "EDYN", the reference number, the lot number and, if available, the serial number and date of manufacture). Read the human-readable interpretation of the bar code symbol printed below for verification.

Training

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System and Dyna (Direct) System.

Copyright and trademarks

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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Dyna Dental Engineering B.V.