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ENGLISH - Istruzioni for use DYNA HELIX® IMPLANTS

Warning : The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Implant Systems. Guidance in the handling of the Dyna Helix® Implant System by an experienced operator is strongly recommended. Dyna Helix® Implant System must only be used by properly trained dentists/doctors and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products.

For more detailed information please refer to the Dyna Implant Manuals as well as Dyna Terms of Guarantee – available on request.

With the publication of this instructions for use all previous versions are no longer valid.

Product description

The Dyna Helix® Implant System consists of titanium cylindrical, screw type implants with ART (Acid etched Rough Titanium) coating and a range of prosthetic abutments. The implants are available in various lengths and diameters. The Dyna Helix® Implants can be used in a one or two stage surgical procedure taken the bone quality and the operation technique in consideration. The system also contains surgical, prosthetic and laboratory components.

Materials

The following materials Dyna uses for its implant products: Titanium (90%Ti-6%Al-4%V).

Content package: See label on packaging.

Indications

The Dyna Helix® Implant System is intended for restoring edentulous or partially edentulous patients with fixed, removable or fixed detachable prostheses. Only use the Dyna Helix® DC implant 6mm for support in combination with several long implants. Dyna Helix® TM implants are specially designed for single stage procedures. The Dyna Helix® Implant System is intended to be used for all types of oral implant supported restorations including single missing teeth. The Dyna Helix® Implant System is fitted with the Dyna Octalock® connection, therefore all prosthetic components are interchangeable with all Dyna Octalock® Implant Systems.

Contraindications

All contraindications associated with elective oral surgery should be taken into consideration. These include: local contraindications (e.g. local infection, inadequate hygiene, unfavourable bite relation, mucosa infections, jaw defects, bruxism, periodontal disease etc), systemic contraindications (e.g. bone vascular disorders, cardiac disease, nephritis/es or nephrosis, rheumatic disease, defective immune response, allergic diseases etc), age related limitations (e.g. children. With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented) and other (e.g. pregnancy, chemotherapy, lack of motivation, diseases treated with large doses of steroids etc.). The use of the Dyna Helix® Implant System is contraindicated, in all situations where it is not indicated.

Surgical method (shortened version)

- ▶ Make sure that the indications for implant placement are met.
- ▶ According to clinical situation, use anaesthetics, make incision and raise the mucosal flap.
- ▶ Under excessive cooling (by sterile saline solution) prepare the implant site using depth markings (6, 8, 10, 11,5, 13, 15, 17 mm). Depths are clearly indicated on (all) drills.
- ▶ Fixed sequence of drilling is obligatory (Helix®: pilot, ø3,2; ø3,6; ø4,2. ø5,0 mm). The recommended speed for the implant site preparation with Dyna drills is max. 800 rpm (decrease the drilling speed when using larger diameter drills: ø3.6-600rpm, ø4.2-500rpm, ø5.0-400rpm). Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the mini-um bone width for ø3.6 implants is 5,6mm).
- ▶ Finalize the shape of the site with the Dyna Helix® Cortical Reamer (Not applicable for the Dyna Helix® DC 6mm, Dyna Helix® TM and Dyna Helix® ST implants). Measure depth according to the bone quality (Type 1= 6mm, Type 2=4-6mm, Type 3=2-4mm, Type 4=0-2mm). For more details see the Dyna Helix® Implant manual.
- ▶ According the bone quality use the Dyna Helix® Taper. In type I bone the Taper should always be used. For more details see the Dyna Helix® Implant manual.
- ▶ Check the site with the Dyna Parallel/depth Instruments/Gauges and perform radiographic examination.
- ▶ Open the implant packaging. The implant is freely seated in the titanium tube (closing screw is included in the lid). Push the Dyna Octa driver for the Torque Wrench or Latch head into the internal Octagon of the implant and remove the implant gently. Insert the implant gently into the preparation; thread it into its final position with a maximum of 30 rpm and an optimal torque range between 25 and 55Ncm, but never exceeding 60Ncm. Tighten the closing screw with 15Ncm or when indicated an abutment with 30Ncm (note: extension abutments with 35Ncm).
- ▶ Prevent (over)loading Dyna Helix® (TM) implants for minimally 8 weeks. When biomechanical criteria for immediate loading are not met this period of time must be extended.
- ▶ Suture the site. Give the patient instructions. If required or necessary make a temporary structure. In order to check the implants' position and control its healing process it is recommended to perform radiographic examination during and/or after the surgery, and later during the healing period. Implants must achieve sufficient primary stability.

Precautions

Proper planning is essential for successful implant treatment.

Beware during the implant planning phase that the marks on the drills for the Helix® implants are situated 0,25 mm higher than the indicated implant length.

Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.

Avoid thermal trauma to the bone; use the largest possible diameter of implants.

Use only in combination with original instruments.

Implant mobility, bone loss, infection may be the symptoms of implant failure.

Helix® Implant should be inserted with the optimal torque ranging between 25-55Ncm.

Mechanical damage of the Helix® Implant can occur above a torque of 70Ncm and of latch-head instruments above a torque of 50Ncm

Avoid overheating during insertion of the implant.

Dyna implants must not be altered in any way.

Rotary instruments must not be used more then 20 times, or after they become damaged or dull.

Despite undergoing special treatments for wear resistance instruments can wear after time. Change the instruments for newer ones if there is any doubt regarding the fit to prevent damaging implants, abutments, screws, etc.

When treating an implant patient, facilitate a disinfected or sterile environment.

Facilitate safe disposal of devices according to safety precautions in the general dental practices or hospitals.

In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

Sterilization

All Dyna implants are sterile and intended for single use only. Under no circumstances sterilize or re-sterilize implants. Do not use implants when the packaging is damaged.

Complications

Apart from the general complications associated with any oral surgery the following complications, related to dental implants, are the most frequent: inflammation, infection, bone loss, swelling, (chronic) pain, paresthesia, patient discomfort, tissue degeneration, bone/implant/restoration fracture, implant mobility, exfoliation and injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Implant Manuals.

Please note:

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any iability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority. When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested.

Handling and Storage

Store in clean, dry, dust-free, dark environment at room temperature. Do not use after expiry date as indicated on the packaging.

Delivery

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

Traceability of serial/lot numbers

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) "EDYN", the reference number, the lot number and, if available, the serial number and expiry date). Read the human-readable interpretation of the bar code symbol printed below for verification.

Training

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System.

Copyright and trademarks

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authori-zation of Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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Copyright en merken: Alle Dyna documenten, zowel geheel als gedeeltelijk, mogen niet gekopieerd, opnieuw gedrukt, bewerkt of gepubliceerd worden zonder schriftelijke toestemming van Dyna Dental Engi- neering BV. Dyna®, Muchor®, Dyna Octalock® Instant Adjusting Bar® en Dyna Helix® zijn geregistreerde merken van Dyna Dental Engineering.