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ENGLISH - Instructions for use DYNA HELIX® IMPLANTS

**Warning**

The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Implant Systems. Dyna Helix® Implant System must only be used by properly trained dentists/dental surgeons and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the Dyna Implant Manuals as well as Dyna Terms of Guarantee – available on request. With the publication of this instructions for use all previous versions are no longer valid.

**Product description**

The Dyna Helix® Implant System consists of titanium cylindrical, screw type implants with ART (Acid etched Rough Titanium) coating. The implants are available in various lengths and diameters. The Dyna Helix® Implants can be used in a one or two stage surgical procedure taken the bone quality and the operation technique in consideration. The system also contains surgical, prosthetic and laboratory components. The Dyna Helix® Implant System is fitted with the Dyna Octalock® connection, therefore all prosthetic components are interchangeable with all Dyna Octalock® Implant Systems.

**Materials**

The following materials Dyna uses for its implant products: Titanium (90%Ti-6%Al-4%V).

**Content package:** See label on packaging.

**Intended purpose**

The Dyna Helix® Implant System is intended for restoring edentulous or partially edentulous patients with fixed, removable or fixed detachable prostheses. The Dyna Helix Implants are intended to be used as anchorage for prosthetic tooth replacement.

**Indications**

The Dyna Helix® Implant System is intended to be used for all types of oral implant supported restorations including single missing teeth. Only use the Dyna Helix® DC implant 6mm for support in combination with several long implants. Dyna Helix® TM implants are specially designed for single stage procedures.

**Contraindications**

All contraindications associated with elective oral surgery should be taken into consideration. In general using implants is contraindicated for all patients who are medically unfit for an oral surgical procedure, patients in whom desired sizes, numbers or positions of implants are not feasible to safely support functional or possibly parafunctional loads and for patients known to be allergic to one or more of the elements in the used material(s). With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented.

**Patient Population**

The patient population consists of adult patients, excluding patients with contraindications associated with elective oral surgery.

**Surgical method (shortened version)**

- Make sure that the indications for implant placement are met.
- According to clinical situation, use anaesthetics, make incision and raise the mucosal flap.
- Under excessive cooling (by sterile saline solution) prepare the implant site using depth markings (6, 8, 10, 11,5, 13, 15mm). Depths are clearly indicated on (all) drills.
- Fixed sequence of drilling is obligatory (ø2,0/2,3, ø3,2, ø3.6, ø4.2, ø5,0 mm). The recommended speed for the implant site preparation with Dyna drills is max. 800 rpm (ø2,0/2,3) (decrease the drilling speed when using larger diameter drills: ø3.2 600-650rpm, ø3.6 400-500rpm, ø4.2 200-350rpm, ø5.0 200rpm). Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the minimum bone width for ø3.6 implants is 5.6mm).
- Finalize the shape of the site with the Dyna Helix® Cortical Reamer (Not applicable when there is lack of high density cortical bone, and for the Dyna Helix® DC 6mm, Dyna Helix® TM and Dyna Helix® ST implants). Measure depth according to the markings. For more details see the Dyna Helix® Implant manual.
- According to the bone quality use the Dyna Helix® Taper. In type I bone the Taper should always be used. For more details see the Dyna Helix® Implant manual.
- Check the site with the Dyna Parallel/depth Instruments/Gauges and perform radiographic examination.
- Open the implant packaging. The implant is freely seated in the titanium tube (closing screw is included in the lid). Push the Dyna Octa driver for the Torque Wrench or Latch head into the internal Octagon of the implant and remove the implant gently. Insert the implant gently into the preparation; thread it into its final position with a maximum of 30 rpm and an optimal torque range between 25 and 55Ncm, but never exceeding 60Ncm. Tighten the closing screw with 15Ncm or when indicated an abutment with 35Ncm.
- Prevent (over)loading Dyna Helix® (TM) implants for minimally 8 weeks. When biomechanical criteria for immediate loading are not met this period of time must be extended.
- Suture the site. Instruct the patient regarding oral hygiene. If required or necessary make a temporary structure. In order to check the implants' position and control its healing process it is recommended to perform radiographic examination during and/or after the surgery, and later during the healing period. Implants must achieve sufficient primary stability.

**Precautions**

- All Dyna implants are sterile and intended for single use only. Under no circumstances sterilize or re-sterilize implants.
- Do not use implants when the packaging is damaged. Contaminated implants must not be used.
- Check the implants and instruments on defects before use.
- Adequate preoperative planning including imaging diagnostics and detailed anamnesis and possible intolerances is essential for successful implant treatment.
- When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested. Secure small parts during operation to protect the patient from aspiration.
- Beware during the implant planning phase that the marks on the drills for the Helix® implants are situated 0,25 mm higher than the indicated implant length.
- Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.
- Avoid thermal trauma to the bone.
- Use only in combination with original instruments.
- Use instruments with shank in combination with contra-angle handpieces for shanks according to ISO 1797 (dental shank).
- Implant mobility, bone loss, infection may be the symptoms of implant failure.
- Helix® Implant should be inserted with the optimal torque ranging between 25-55Ncm.
- Mechanical damage of the Helix® Implant can occur above a torque of 70Ncm and of latch-head instruments above a torque of 50Ncm
- Avoid overheating during insertion of the implant.
- Dyna implants must not be altered in any way.
- Check the instruments on defects before use. Instruments have a limited number of reprocessing cycles. Discard defective instruments (faded depth marking, blunt/chipped blades, bent/fractured working parts, corroded surfaces).
- When treating an implant patient, facilitate a disinfected or sterile environment.
- Facilitate safe disposal of devices according to safety precautions in the general dental practices or hospitals. Decontaminate contaminated products/instruments if applicable.
- In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

**MRI safety information**

Please note that these products have not been evaluated for safety and compatibility in the MR environment. These products have not been tested for heating or migration in the MR environment. For additional information on Magnetic Resonance Imaging, please consult the MRI statement available at www.dynadental.com/IFU/MRI.

**Complications**

Apart from the general complications associated with any oral surgery the following complications, related to dental implants, are the most frequent: inflammation, infection, bone loss, swelling, (chronic) pain, paresthesia, patient discomfort, tissue degeneration, bone/implant/ restoration fracture, implant mobility, exfoliation and injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Helix® Implant manual.

**Please note:**

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

**Handling and Storage**

Fragile; handle with care. Store in clean, dry, dust-free, dark environment at room temperature. Do not use after expiry date as indicated on the packaging.

**Delivery**

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

**Traceability of serial/lot numbers**

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) "EDYN", the reference number, the lot number and, if available, the serial number and expiry date). Read the human-readable interpretation of the bar code symbol printed below for verification.

**Training**

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System.

**Copyright and trademarks**

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Dyna Octalock® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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**Copyright en merken:** Alle Dyna documenten, zowel geheel als gedeeltelijk, mogen niet gekopieerd, opnieuw gedrukt, bewerkt of gepubliceerd worden zonder schriftelijke toestemming van Dyna Dental Engineering BV. Dyna®, Dyna Octalock® en Dyna Helix® zijn geregistreerde merken van Dyna Dental Engineering.