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ENGLISH - Instructions for use cleaning_sterilisation_reprocessing

Scope of the IFU

These instructions for use contain the instructions for initial cleaning, disinfection and sterilisation of the Abutments, Attachments, Instruments and Implant Products (hereafter: components) and for reprocessing of reusable products. For instructions on using the medical devices manufactured by Dyna Dental Engineering B.V., please refer to the corresponding instructions for use.

Read these instructions for use completely and carefully. Keep the instructions for use for future reference. Injuries to the patient or damage to the medical device cannot be ruled out if the instructions for use are not observed.

⚠ The washer-disinfectors and small steam sterilizers used must be regularly maintained and checked. All validated parameters must be complied with at all times.

⚠ Magnets must only be cleaned and disinfected. These products must not be sterilized, as they loose their intended purpose due to the high temperature effect.

Cleaning, disinfection and sterilisation

Effective cleaning and disinfection are the necessary prerequisites for sterilisation for all components.

⚠ The components are supplied non-sterile and must be cleaned, disinfected and sterilized by the user before the first use.

The user is responsible for ensuring that the cleaning, disinfection and sterilisation actually carried out, using equipment, materials and personnel used in these procedures, achieves the desired results. This usually requires qualifications and routine monitoring of the materials and equipment used.

Preparation of cleaning

The components can be cleaned without further preparation and disassembly. All products must be treated individually before cleaning and disinfection.

Manual Pre-cleaning

It is recommended to process the instruments within one hour of use at very latest.

⚠ Pre-cleaning does not replace automated cleaning and disinfection and serves, among other things, to protect personnel.

⚠ Cleaner for instruments must not contain any aldehydes or alcohols. The concentrations and minimum take effect times specified by the manufacturer and the instructions for rinsing must be complied with.

- Place components in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Comet DC Evo cleaning and disinfecting agent) to prevent drying of residues on the products (protein fixation) and to facilitate the cleaning of the products.
- Remove components from cleaning/disinfection tank and rinse components thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the machine.

Manual pre-cleaning is directly followed by mechanical cleaning and disinfection.

Mechanical cleaning and disinfection

⚠ Only use a standard-compliant washer-disinfecter according to EN ISO 15883-1 with subsequent product drying.

⚠ Use enzymatic cleaning concentrate for Thermo washer disinfectors e.g. Euroclean 120 and Rinse aid concentrate for Thermo washer disinfectors e.g. Eurobright 360

- Place the components in a suitable bur block or small parts holder.
- Place the bur block/small parts holder in the washer/disinfecter in such a way that the products are directly hit by the spray jet.
- Start the Vario TD programme or universal programme (for parameters of program see Table 1) including thermal disinfection. Thermal disinfection takes play allowing for the A₀ value ≥ 3000.
- On completion of the cycle remove components from the washer/disinfecter.
- Visual examination to ensure that the products are clean and undamaged.

⚠ If after mechanical cleaning and disinfection there are still visible residues of contamination, repeat the cleaning and disinfection process until no visible contamination is left.

Repeat steps 1-5 if the products aren't clean until no visible contamination is left.

The components must be packed immediatly after cleaning and disinfection.

Table 1: Cleaning and disinfection parameters

| Pre-cleaning | Cleaning | Disinfection | | A ₀ value |
|--------------|----------|---------------------|-------------|----------------------|
| | | Duration | Temperature | |
| < 45 °C | ≥ 55 °C | ≥ 300 s (5 minutes) | ≥ 90 °C | ≥ 3.000 |

Packaging and Sterilisation

⚠ All components must not be sterilized in the original packaging. Use an appropriate seal-method to pack the instruments in sterilisation foil.

⚠ Magnets must not be sterilized.

- Make sure that the packaging is suitable for the components and the chosen method of sterilisation. The packaging must be large enough to ensure that there is no pressure on the seal. Use an appropriate method to pack the components. Products with limited use are to be marked accordingly.
- Place the components onto a suitable sterilisation tray. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several products.
- Start the sterilisation program (for parameters see Table 2).

Table 2: Sterilisation parameters

| Temperature | Pressure | Holding time | Drying time |
|-------------|------------|--------------|-------------|
| 134 °C ± 1% | 3 bar ± 5% | ≥ 3 minutes | ≥ 5 minutes |

Storage

After sterilisation the packed sterile goods must be protected from dust, moisture and recontamination during transport and storage and labeled with the sterilisation date, contents and release after successful sterilisation.

Disposal

Disposal of the components must be carried out in accordance with local regulations and environmental requirements, considering the respective degree of contamination. Cleaning as well as disinfection before disposal may be necessary.

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