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ENGLISH - Instructions for use DYNABUTMENTS

Warning

The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Implant Systems. Dyna Helix® Implant Systems must only be used by properly trained dentists/dental surgeons and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the Dyna Implant Manuals as well as Dyna Terms of Guarantee – available on request. With the publication of this instructions for use all previous versions are no longer valid.

Product description

The Dyna (Octalock®) abutments are a part of the Dyna (Helix®, Octalock®) Implant System, and used in the prosthetic part of implant rehabilitation. All abutments (except for the magnetic-Dyna alloy; and PEEK abutments/caps) are made out of titanium. They are meant to cover all clinical and prosthetic situations. Three different abutment lines are supplied by Dyna:

- Dyna abutments - to be used for the Dyna push-in implants (no octa).
- Dyna Octalock® abutments - to be used as well for the Dyna Octalock® push-in implants as the Dyna Helix® Implants.
- Dyna competitor screws - to be used for CadCam constructions on competitor implants.

The following materials Dyna uses for its abutments: Titanium (90%Ti-6%Al-4%V), Precious alloy (60%Pd-36,8%Co-1%Pt)-2,2%Other elements (Ga,Cr, etc.), PEEK (Polyetheretherketone). Products must not be used for patients known to be allergic to one or several of the elements contained in the abutment materials. In case of doubt perform preliminary allergological testing.

Content package

See label on packaging.

Intended purpose

Dyna abutments in combination with endosseous implants are used as the foundation for anchoring tooth replacements in either jaw and are intended to restore function and aesthetics in the oral cavity.

Indications

- Dyna ball, instant adjusting bar or magnet abutments are intended to be used for removable (overdentures) prostheses on edentulous or partially edentulous patients.
- Dyna temporary, straight, angled, tibase and anatomic abutments are intended to be used for fixed constructions restoring single or multiple missing teeth on partially edentulous patients.
- Dyna extension and multi unit abutments are intended to be used to elevate seating platform of restorations.
- Multi Unit abutments in particular are intended to be used for fixed detachable prostheses for multiple unit (minimal 4) screw-retained restoration and to overcome angulations up to 30°.

Contraindications

All contraindications associated with elective oral surgery and prosthetic rehabilitation on implants should be noticed. In general using abutments is contraindicated for all patients who are medically unfit for an oral surgical procedure, patients in whom desired sizes, numbers or positions of implants are not feasible to safely support functional or possibly parafunctional loads and for patients known to be allergic to one or more of the elements in the used material(s). With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented).

The use of Dyna abutments for fixed constructions (temporary, straight, angled, tibase and anatomic abutments) is contraindicated:

- for individual tooth restorations with free end saddle.
- or restorations with angulation correction of more than 30° to the implant axis.

The use of Dyna abutments for removable constructions (ball and magnet) is additionally contraindicated in situations where:

- proper articulation cannot be restored
- proper alignment of the male and female parts cannot be guaranteed

The use of the Dyna Instant Adjusting Bar is additionally contraindicated for:

- implants placed less than 16 mm or more than 26mm from each other
- implants with an angulation of more than 9° in a non parallel position
- in situations when occlusal forces transferred to the bar are not favourable
- unfavourable bucco/lingual supporting implants inclination

Patient Population: The patient population consists of adult patients, excluding patients with contraindications associated with elective oral surgery and prosthetic rehabilitation.

Making impression

For Dyna Octalock® System (fixed and bar constructions) and for Dyna Implant System (bar constructions) use open-tray technique to transfer the position of the implant to the working model in the lab. For all other constructions use either special transfer elements or conventional approach to make an impression. See the Dyna Prosthetics manuals for detailed information.

Abutment choice

The diameter, angulation of the implant, the gingival height, hygienic requirements and patients dexterity must be taken into consideration when selecting the abutment height. If an improper abutment is used the final prosthetic result may be challenged or it may lead to the implant failure.

Antirotation

Dyna Octalock® System has a special anti-rotation system preventing unscrewing of the abutment from the implant. All-in-one abutments from Dyna Implant System should be given special care to prevent unscrewing (see also Dyna Implant Manual). Always use a torque instrument and tighten the abutments with 35Ncm torque force (note: cover and healing screws and the Multi-Unit Prosthetic Screw Hexagon/ CADCAM Bar screw Dyna Multi-Unit with 15Ncm).

Procedure

See detailed information regarding all meso- and supra structure possibilities in the Dyna manuals. For detailed information about tooling and torque of competitor abutments or screws see their manuals.

Precautions

- Check the abutments on defects before use.
- Adequate preoperative planning including imaging diagnostics and detailed anamnesis and possible intolerances is essential for successful implant treatment.
- To prevent screw loosening or fracture, retighten a screw/abutment on implant level with 35Ncm before the final prosthetic treatment.
- When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested. Secure small parts during operation to protect the patient from aspiration
- Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.
- Use only in combination with original instruments.
- Implant mobility, bone loss, infection may be the symptoms of implant failure.
- When treating an implant patient, facilitate a disinfected or sterile environment.
- Facilitate safe disposal of devices according to safety precautions in the general dental practice or hospitals. Decontaminate contaminated products/instruments if applicable.
- In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer

Sterilization

All Dyna abutments are supplied not sterile in a blister packaging, and intended for single use only. They must be disinfected and sterilized before clinical use as stated in the IFU cleaning_sterilisation_reprocessing available at www.dynadental.com/IFU/. Do not use abutments when the packaging is damaged.

Handling and Storage

Store in ambient conditions.

MRI safety information

Please note that these products have not been evaluated for safety and compatibility in the MR environment. These products have not been tested for heating or migration in the MR environment. For additional information on Magnetic Resonance Imaging, please consult the MRI statement available at www.dynadental.com/IFU/MRI.

Complications

Possible complications, related to prosthetic rehabilitation with dental implants, include: inflammation, bone loss, (chronic) pain, patient discomfort, tissue degeneration, bone/implant/restoration fracture, implant/suprastructure mobility. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Helix® Implant manual.

Please note

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

Delivery: Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

Traceability of serial/lot numbers: It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) «EDYN», the reference number, the lot number and, if available, the serial number and date of manufacture). Read the human-readable interpretation of the bar code symbol printed below for verification.

Training: Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System.

Copyright and trademarks: All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Dyna Octalock® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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Copyright en merken: Alle Dyna documenten, zowel geheel als gedeeltelijk, mogen niet gekopieerd, opnieuw gedrukt, bewerkt of gepubliceerd worden zonder schriftelijke toestemming van Dyna Dental Engineering BV. Dyna®, Dyna Octalock® en Dyna Helix® zijn geregistreerde merken van Dyna Dental Engineering.