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ENGLISH - Instructions for use DYNAATTACHMENTS

Warning

The descriptions given in this enclosure are insufficient to allow immediate use of the Dyna (Direct) System and Dyna Helix® Implant System. Both systems must only be used by properly trained dentists/dental surgeons and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the the Dyna Magnet and Dyna Helix® Implant Manual as well as Dyna Terms of Guarante – available on request. With the publication of this instructions for use all previous versions are no longer valid.

Product description

Dyna attachments include the magnetic attachment system – magnetically retained prostheses consisting of magnets and a pre-fabricated, ferromagnetic alloy keeper to be used directly on endodontically treated roots.

Materials

Dyna uses for its magnets: case of Surgical Steel (316L) and magnet 26%Nd-70%Fe-1%B (Zn coated), Dyna uses for its keepers: Precious alloy (60%Pd-36,8%Co-1%Pt)-2,2%Other elements (Ga,Cr, etc.),

Content package

See label on packaging.

Intended purpose

Dyna magnetic systems are intended for increasing retention of (partial) overdentures, supported on implants and natural roots. Magnets are specially indicated for overdentures in patients with manual dexterity deficiency and/or severely reduced alveolar bone.

Indications

Dyna magnetic systems are intended for increasing retention and stabilisation of (partial) overdentures, supported on implants and natural roots (also for elements with unfavourable prognosis). Magnets are specially indicated for overdentures in patients with manual dexterity deficiency and/or severely reduced alveolar bone. Prefabricated keepers (Dyna Direct System) are to be used for elements with a favourable inclination and in case the preparation can be placed above the gingiva level.

Contraindications

All contraindications associated with elective oral surgery and prosthetic rehabilitation on implants and endodontically compromised teeth should be noticed. In general using attachments is contraindicated for all patients who are medically unfit for an oral surgical procedure, patients in whom desired sizes, numbers or positions of implants or endodontically compromised teeth are not feasible to safely support functional or possibly parafunctional loads and for patients known to be allergic to one or more of the elements in the used material(s). With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented).

The use of the Dyna magnetic system is additionally contraindicated in situations where proper articulation can not be restored and proper alignment of the male and female parts can not be guaranteed.

Patient Population

The patient population consists of adult patients, excluding patients with contraindications associated with elective oral surgery and prosthetic rehabilitation, also on endodontically compromised teeth.

Procedure

Each product has its own special clinical and laboratory procedure. See the Dyna Magnet and Dyna Helix® Implant Manual for detailed information.

Clinical application

Whenever working with overdenture attachments take the following into consideration:

- ▶ allowing the patient 2-4 weeks of wear-in time
- ▶ type of implant/root-overdenture connection (resilient vs. non-resilient)
- ▶ alignment of the male and female parts

For details see the Dyna Magnet and Dyna Helix® Implant Manual.

Precautions

- Check the attachments on defects before use. Do not grind or sandblast.
- Improper technique can contribute to implant/prosthesis failure and/or bone loss.
- Check the overdenture/attachments regularly (every 6 months).
- When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested. Secure small parts during operation to protect the patient from aspiration
- Use only in combination with original instruments.
- Denture cleaning agents should be used with care. Overnight cleaning solutions should be avoided.
- Implant mobility, bone loss, infection may be the symptoms of implant failure.
- Magnets may not be heated above 120°C / 248°F and only be used in combination with Dyna EFM alloy, Dyna Keepers and Dyna Medical(magnet) abutments.
- Denture without balanced occlusion and articulation will contribute to failure.
- When treating an implant patient, facilitate a disinfected or sterile environment.
- Facilitate safe disposal of devices according to safety precautions in the general dental practice or hospitals. Decontaminate contaminated products/instruments if applicable.
- In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

Sterilization

All Dyna attachments are supplied not sterile in a blister packaging, and intended for single use only. They must be disinfected and sterilized before clinical use as stated in the IFU cleaning_sterilisation_reprocessing available at www.dynadental.com/IFU/. Magnets must only be cleaned and disinfected. These products must not be sterilized, as they loose their intended purpose due to the high temperature effect. Do not use attachments when the packaging is damaged.

Handling and Storage

Store in ambient conditions.

Magnetic fields

The research on influence of the magnetic fields on mammalian tissues shows that using magnets for dental purposes poses no threat to the patients. The influence of the magnets on pacemakers devices is limited and when magnets are used according to their indications negligible.

MRI safety information

Please note that these products have not been evaluated for safety and compatibility in the MR environment. These products have not been tested for heating or migration in the MR environment. Magnets may influence the image of MRI. Magnets should therefore always be removed before the MRI scan. For additional information on Magnetic Resonance Imaging, please consult the MRI statement available at www.dynadental.com/IFU/MRI.

Complications

Possible complications, related to prosthetic rehabilitation with attachments on natural roots include:

inflammation, (chronic) pain, patient discomfort, root fracture, increase in root mobility, perforation of the root. In case of complications follow the specific for the particular situation, appropriate course of actions - for more details see the Dyna Magnet Manual.

Please note:

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

Delivery

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

Traceability of serial/lot numbers

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) “EDYN”, the reference number, the lot number and, if available, the serial number and date of manufacture). Read the human-readable interpretation of the bar code symbol printed below for verification.

Training

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System and Dyna (Direct) System.

Copyright and trademarks

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Dyna Octalock® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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Dyna Dental Engineering B.V.