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**ENGLISH - Instructions for use** **INSTRUMENTOS DYNÁ**

**Warning**

The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Systems. Dyna Helix® Implant System must only be used by properly trained dentists/dental surgeons and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the Dyna Magnet and Dyna Helix® Implant Manual as well as Dyna Terms of Guarantee – available on request. With the publication of this instructions for use all previous versions are no longer valid.

**Product description**

The Dyna instruments group contains:

- a range (pilot, spiral drills and auxiliary) of specially designed drills for the Dyna Helix® implants, available in various lengths and diameters, as well as drills for the Dyna Direct System. All drills are externally irrigated.
- cortical reamers and tapers for the Dyna Helix® implants, used for finalizing the shape of the bone site before introducing the implant.
- acta drivers for the Dyna Helix® implants, used for seating the implants and thread these into its final position.
- sulcus reamer for the Dyna Helix® implants, used for shaping the supra-implant part of the bone crest, before abutment tightening.

All instruments are used in combination with surgical handpieces.

**Materials**

Dyna instruments are made out of Stainless Steel (DIN 1.4112, 1.4197).

**Content package**

See label on packaging.

**Intended purpose**

The Dyna instruments are intended to be used for preparation of the bone site and sealing of the implants during intra-oral implant treatment for Dyna Implant Systems and/or oral rehabilitation in combination with Dyna products and the natural roots for the Dyna Direct System.

**Indications**

The Dyna instruments are intended to be used for preparation of the bone site and sealing of the implants during intra-oral implant treatment for Dyna Implant Systems and/or oral rehabilitation in combination with Dyna products and the natural roots for the Dyna Direct System.

- Pilot and Spiral drills are used to prepare the bone site for the Dyna Helix Implants.
- The Dyna Direct Spiral and Seatdrill are used to drill the necessary preparation for the Dyna Direct keeper.
- Cortical reamers and tapers for the Dyna Helix® implants are intended to be used for finalizing the shape of the bone site before especially in high density bone.
- Octa drivers for the Dyna Helix® implants are intended to be used for seating the implants and thread these into its final position.
- The sulcus reamer for the Dyna Helix® implants is intended to be used for shaping the supra-implant part of the bone crest, before abutment tightening.

All instruments are intended to be used in combination with surgical handpieces.

**Contraindications**

All contraindications associated with elective oral surgery should be noticed. In general using abutments is contraindicated for all patients who are medically unfit for an oral surgical procedure, patients in whom desired sizes, numbers or positions of implants are not feasible to safely support functional or possibly parafunctional loads and for patients known to be allergic to one or more of the elements in the used material(s). With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented). Cortical reamers should not be used when there is lack of high density cortical bone, and for the Dyna Helix DC Ø 4.2 L 6mm, the Dyna Helix ST and TM implants. Tapers should not be used for the Dyna Helix DC Ø 4.2 L 6mm implant.

**Patient Population**

The patient population consists of adult patients, excluding patients with contraindications associated with elective oral surgery.

**Preparation of the implant bed**

- Make sure the bone quality and quantity is adequate and that the indications for implant placement are met.
- According to clinical situation, use anaesthetics, make incision and raise the mucosal flap.

- Under excessive cooling (by sterile saline solution) prepare the implant site using depth markings (6, 8, 10, 11.5, 13, 15mm). Depths are clearly indicated on (all) drills.

- Fixed sequence of drilling is obligatory (ø2.0/2.3, ø3.2, ø3.6, ø4.2, ø5.0 mm). The recommended speed for the implant site preparation with Dyna drills is max. 800 rpm (ø2.0/2.3) (decrease the drilling speed when using larger diameter drills: ø3.2 600-650rpm, ø3.6 400-500rpm, ø4.2 200-350rpm, ø5.0 200rpm. Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the minimum bone width for ø3.6 implants is 5.6mm).

- Finalize the shape of the site with the Dyna Helix® Cortical Reamer (Not applicable when there is lack of high density cortical bone, and for the Dyna Helix® DC 6mm, Dyna Helix® TM and Dyna Helix® ST implants). Measure depth according to the markings. For more details see the Dyna Helix® Implant manual.

- According to the bone quality use the Dyna Helix® Taper. In type I bone the Taper should always be used. For more details see Dyna Helix® Implant manual.

- Check the site with the Dyna Parallel/depth Instruments/Gauges and perform radiographic examination.

**Preparation of the root for Dyna Direct System**

- On the chosen teeth perform an endodontic treatment in an usual manner. Decoronate the abutment tooth up till 1mm subgingival. Prepare the final surface more or less parallel to the alveolar crest but with a slight lingual or palatinal inclination

- Drill a shaft in the root using the Dyna spiral drill and prepare the seat using the Dyna seat drill.

**Precautions**

- Adequate preoperative planning including imaging diagnostics and detailed anamnesis and possible intolerances is essential for successful implant treatment.

- When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested. Secure small parts during operation to protect the patient from aspiration.

- Beware during the implant planning phase that the marks on the drills for the Helix® implants are 0,25 mm higher as the indicated implant length.

- Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.
- Use only in combination with original instruments and Dyna implants.
- Use instruments with shank in combination with contra-angle handpieces for shanks according to ISO 1797 (dental shank).
- Check the instruments on defects before use. Instruments have a limited number of reprocessing cycles. Discard defective instruments (faded depth marking, blunt/chipped blades, bent/fractured working parts, corroded surfaces).
- Mechanical damage of latch-head instruments can occur above a torque of 50Ncm. Then use a T.W. instrument/driver in combination with a hand-used torque wrench.
- Avoid thermal trauma to the bone – sufficient cooling.
- When treating an implant patient, facilitate a disinfected or sterile environment.
- Facilitate safe disposal of devices according to safety precautions in the general dental practice or hospitals. Decontaminate contaminated products/instruments if applicable.
- In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

**Sterilization**

All instruments are supplied not sterile in a blister packaging. They must be disinfected and sterilized before clinical use as stated in the IFU cleaning\_sterilisation\_reprocessing available at www.dynadental.com/IFU/. Do not use instruments when the packaging is damaged.

**Handling and Storage**

Store in ambient conditions.

**Complications**

Possible complications, related to intra-oral bone site preparation with implant drills include: infection, bone/tissue/nerve damage, thermal trauma, excessive bleeding, sinus perforation, patient discomfort, tissue degeneration, drill fracture, injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery. For more details see the Dyna Helix® Implant manual.

**Please note:**

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

**Delivery**

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

**Traceability of serial/lot numbers**

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) “EDYN”, the reference number, the lot number and, if available, the serial number and date of manufacture). Read the human-readable interpretation of the bar code symbol printed below for verification.

**Training**

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System.

**Copyright and trademarks**

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Dyna Octalock® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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