ENGLISH - Instructions for use DYNA instruments (class IIa)

**Warning**

The descriptions given in this enclosure are insufficient to allow immediate use of all Dyna Systems. Guidance in the handling of the Dyna (Direct) System, Dyna Helix Implant System by an experienced operator is strongly recommended. Dyna Helix® Implant Systems must only be used by properly trained dentists/doctors and in combination with original components. In case of multiple use the following dangers could occur: cross infection and/or damaged products. This can result in loss of function and wrong identification of products. For more detailed information please refer to the Dyna Manuals as well as Dyna Terms of Guarantee – available on request. With the publication of this instructions for use all previous versions are no longer valid.

**Product description**

Dyna instruments (e.g. latch-head) are class IIa medical products. They contain range (pilot, spiral drills and auxiliary) of specially designed drills for the Dyna Helix® implants, available in various lengths and diameters-and Dyna Direct System drills. All drills are externally irrigated. As a rule Dyna drills are used in combination with surgical handpieces. Other instruments are tapers,-reamers (e.g. cortical reamer) and drivers.

**Materials:** Dyna instruments are made out of Stainless Steel (DIN 1.4034, 1.4108, 1.4112, 1.4197).

**Content package:** See label on packaging.

**Indications**

The Dyna instruments are intended to be used for preparation of the bone site during intra-oral implant treatment for Dyna Implant Systems and/or oral rehabilitation in combination with Dyna products-and the natural roots for the Dyna Direct System.

**Contraindications**

All contraindications associated with elective oral surgery should be noticed. These include: local (e.g. local infection, inadequate hygiene, unfavorable bite relation, mucosa infections, jaw defects etc), systemic (e.g. bone vascular disorders, cardiac disease, nephritisos or nephroses, rheumatic disease, defective immune response, allergic diseases etc), age related limitations (e.g. children. With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth phase has been properly documented) and other (e.g. pregnancy, chemotherapy, lack of motivation, diseases treated with large doses of steroids etc.). The use of all Dyna instruments is contraindicated in all situations where it is not indicated.

**Preparation of the implant bed**

► Before preparation of the implant bed make sure the bone quality and quantity is adequate and that the indications for implant placement are met.

► Verify the depth of the preparation both pre- and intra-surgically. Performing radiographic examination, using the depth gauges and depth markings is essential. Depth markings for the 6, 8, 10, 11,5, 13 and 15mm implants are clearly indicated on (all) drills.

► Drills must not be damaged, blunt or used more than 20 times.

► Under excessive cooling (sterile saline solution) prepare the implant site; fixed sequence of drilling is obligatory (pilot, ø3.2, ø3.6, ø4.2, ø5.0 mm) the recommended speed for the implant site preparation with Dyna drills is max. 800 rpm (pilot) (decrease the drilling speed when using larger diameter drills: ø3.2 600-650rpm, ø3.6 400-500rpm, ø4.2 200-350rpm, ø5.0-<200rpm

► Finalize the shape of the site with the Dyna Helix® Cortical Reamer (This does not apply to the Dyna Helix® DC 6mm, Dyna Helix® TM and Dyna Helix® ST implants). Depth according to the bone quality. For details see Dyna Helix® Implant manual.

► According to the bone quality use the Dyna Helix® Taper. In type I bone the Taper always should be used. For details see Dyna Helix® Implant manual.

► Placement of Dyna implants requires sufficient bone volume. Implants must always be surrounded by at least 1mm of healthy bone (e.g. the minimum dentoalveolar width for ø3.6 implants is 5.6mm).

► In order to check the preparation position and angulation use the Dyna Parallel/Depth instruments and gauges and perform radiographic examination.

**Preparation of the root for Dyna Direct System**

► On the chosen teeth perform an endodontic treatement in an usual manner. Decoronate the abutment tooth up till 1mm subgingival. Prepare the final surface more or less parallel to the alveolar crest but with a slight lingual or palatal inclination

► Drill a shaft in the root using the Dyna spiral drill and prepare the seat using the Dyna seat drill.

**Precautions**

Proper planning is essential for successful implant treatment.

Beware during implant planning that the marks on the drills for the Helix® implants are 0,25 mm higher as the indicated implant length.

Improper technique can contribute to implant failure and/or bone loss; hard tissues must be treated with care.

Use only in combination with original instruments and Dyna implants.

In spite of special treatments for wear resistance instruments can wear after time. Change the instruments for new when there is any doubt regarding the fit to prevent damaging implants, abutments, screws, etc Mechanical damage of latch-head instruments can occur above a torque of 50Ncm. Then use a T.W. instrument/driver in combination with a hand-used torque wrench.

Avoid thermal trauma to the bone – sufficient cooling.

Facilitate safe disposal of devices according to safety precautions in the general dental practice or hospitals.

In case of any doubts concerning the use of Dyna products contact Dyna Dental Engineering BV or your local dealer.

**Sterilization**

All Dyna instruments are supplied non-sterile and must be sterilized or disinfected in an appropriate manner before clinical use. After each use clean the drills in the following order: running tap water and soft brush, disinfection bath (preferably in ultrasonic washers) and sterilization in an appropriate way.

**Complications**

Possible complications, related to intra-oral bone site preparation with implant drills include: infection, bone/tissue damage, thermal trauma, patient discomfort, tissue degeneration, drill fracture, injury to adjacent anatomic structures. In case of complications follow the appropriate course of actions generally applied in oral surgery -for more details see the Dyna Implant Manual.

**Please note:**

The user of Dyna products is obliged to determine whether or not any products are suitable for a particular clinical situation. The user of Dyna products is obliged to document the products used for each patient in an appropriate manner. Dyna Dental Engineering BV disclaims any liability, express or implied and shall not be responsible for any damages arising from or in connection with any errors in professional judgment or practice in the use or installation of Dyna products. It is the users duty to study the latest developments in dental implantology as well as Dyna Implant Systems and its applications, and to brief the patient about any warnings, precautions, contra-indications, measures to be taken and limitations of use regarding the device. The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority. When using a Dyna product intra-orally take proper care to prevent them from being inhaled or ingested.

**Handling and Storage**

Store in clean, dry, dust-free, dark room at room temperature.

**Delivery**

Federal law restricts these devices to be sold by or on the behest of a dentist or a physician.

**Traceability of serial/lot numbers**

It is the end users responsibility by law to record the serial and/or lot numbers of all products for traceability purposes. The Health Industry Bar Code (HIBC) on the label contains information for traceability purposes (the Labeler Identification Code (LIC) “EDYN”, the reference number, the lot number and, if available, the serial number and date of manufacture). Read the human-readable interpretation of the bar code symbol printed below for verification.

**Training**

Dyna Dental Engineering BV arranges regular training courses for beginning and advanced implantologists. The courses are obligatory and are meant to provide the Dyna user with practical and theoretical expertise concerning the use of the Dyna Implant System.

**Copyright and trademarks**

All Dyna documents may not be copied, reprinted or published in whole or in part without the written authorization of Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® and Dyna Helix® are registered trademarks of Dyna Dental Engineering BV.

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**Copyright en merken:** Alle Dyna documenten, zowel geheel als gedeeltelijk, mogen niet gekopieerd, opnieuw gedrukt, bewerkt of gepubliceerd worden zonder schriftelijke toestemming van Dyna Dental Engineering BV. Dyna®, Muchor®, Dyna Octalock®, Instant Adjusting Bar® en Dyna Helix® zijn geregistreerde merken van Dyna Dental Engineering.
Dyna Dental Engineering B.V. Vang 9